

The World Leading Provider of High Pressure Equipment for Research and Industry since 1945!

100 & 300 ml

Bolted Closure Stirred Laboratory Reactor



At a Glance


Volume: 100 ml & 300 ml

Material of

Construction: 316 Stainless Steel & Hastelloy® C-276

Design Pressure: 5,500 psi @ 650° F
(379 bar @ 343°C)

Applications: The bolted closure stirred laboratory reactor is a versatile high pressure and high temperature, stirred laboratory reactor. It is used for chemical synthesis of corrosive, hazardous and very reactive chemicals / petrochemicals as well as materials research.



Principle of Operation

The Parker Autoclave Engineers' Bolted Closure Reactor has been designed for reliable high pressure operation. The seal is a metal gasket machined from the same material as the vessel. Many combinations of standard components are available. The cover of the unit remains fixed in the stand to permit opening of the vessel without disassembling any process connections. The body is easily removed and drops away from the cover.

Features

- Versatile product configuration
- Operating pressures as high as 4740 psi @ 650° F (327 bar @ 343° C)
- Open vessel and remove body without disassembling pressure connections
- Available worldwide to meet codes such as ASME, CE, and CRN.

General Specifications

Design Pressure

5,500 psi @ 650° F* (379 Bar @ 343° C)*

Minimum Design Metal Temperature (MDMT)

-20° F @ 5,500 psi (-29° C @ 379 Bar)

Maximum Operating Pressure (MOP)

Varies based on gauge, transducer, and rupture disk selection. Refer to Ordering Guide for Details.

Critical Dimensions:	100 ml	300 ml
Inside Diameter:	1.81" (46mm)	1.81" (46 mm)
Straight Wall:	2.75" (70 mm)	6.69" (170 mm)
Approximate Dimensions:	Short Bench Top	Tall Bench Top
Overall Height**	34.8" (883 mm)	39.4" (1002 mm)
Width:	20.5" (522 mm)	20.5" (522 mm)
Depth:	26.2" (665 mm)	26.2" (665 mm)

* 650° F (343° C) rating is vessel mean wall temperature. Actual Process temperature will be lower.

** Overall height based on belt driven units. For actuals see standard drawings.



300 ml Bolted Closure Reactor Internals

Connection Schedule

All of the connections shown will be provided. For any accessories not ordered, the corresponding connection will be plugged.

Opening Label	Purpose	Opening or feature description on underside of cover	Opening or feature description on topside of cover	Entry Point	Smallest diameter orifice (nominal) in flow path
A	Charging Port	0.161" port	3/8" O.D. Tube	Cover Top	0.161"
B	Gas Inlet	Branched into A	1/8" O.D. Tube	Cover Side	0.062"
C	Blow Pipe or Sparge Tube†	1/8" O.D. tube	1/8" O.D. Tube	Cover Side	0.062" 0.031"†
D & H	Cooling Coil	1/8" O.D. tube	1/4" O.D. Tube	Cover Side	0.062"
E	Vent and Pressure Indication	Branched into F	1/8" O.D. Tube	Cover Side	0.062"
F	Safety Head	0.161" port	1/8" FNPT	Cover Top	0.161"
G	Thermowell‡	1/8" O.D. tube	3/32"Port‡	Cover Top	N/A
J	Liquid Sample or Sparge Tube†	1/8" O.D. tube	1/8" O.D. Tube	Cover Side	0.062" 0.031"†
K	MagneDrive® Agitator	1/2" O.D. Mixing Shaft	MagneDrive®	Cover Top	N/A

† The tube that forms the sparge tube is 1/8" O.D. and 1/16" I.D. with a plug in the end. Nine .031" diameter holes are drilled in the sparge ring to bubble gas into the reactor.

‡ The tube that forms the thermowell is 1/8" O.D. and 1/16" I.D. with a plug in the end. A 3/32" port is drilled in the cover to guide the thermocouple to the opening in the thermowell.

Technical Specifications

Parker Autoclave Engineers provides a variety of optional accessories to custom configure each reactor. See the Bolted Closure Stirred Reactor Ordering Guide to configure a reactor for a specific application.

Seal Materials: Metal Gasket (vessel material), Buna-N, Ethylene-Propylene, PTFE, Viton®, Silicone or Kalrez® O-rings

Approvals: Optional ASME code stamp, Canadian Registration or CE Mark

Stand: Short Bench Top or Tall Bench Top

Body Lift: None or Manual Jack

Agitator: 3300 RPM rated MagneDrive® MAG075-01 Series with 7 in-lb (0.79 N-m) static torque. 3300 RPM rated MagneDrive® MAG075-02 Series with 16 in-lb (1.8 N-m) static torque, carbon/graphite bearings or Fluoropolymer with graphite fiber.

Motors: 1/2 HP (0.37 kW) General Purpose DC with either: 90 V Armature (120 V unit), or 180 V Armature (240 V unit) CE Mark. 1/2 HP (0.37 kW) Explosion-Proof DC with either: 90 V Armature (120 V unit), or 180 V Armature (240 V unit). Air Motor with manual or electronic speed adjustment (Supply with 35 SCFM of 40 psi compressed air minimum).

Impeller Styles: AE Dispersimax, Straight Turbine, Axial Flow-Up, or Axial Flow-Down; All 7/8 inch (22.2 mm) diameter.

Baffle: One (1) Single Blade Baffle attached to the top cover is included.

Speed Sensor: General Purpose or Intrinsically-Safe Magnetic Sensor (Barrier Required).

Heating: Furnaces: 120 VAC, Single Phase or 240 VAC, Single Phase; 1,200 Watt. Jacket: Removable, Spiral Baffle with O-Ring Seals.

Internal Accessories Available:

Liquid sample tube w or w/o valve (1/8" O.D. tube)
Blow pipe w or w/o valve (1/8" O.D. tube)
Sparge tube w or w/o valve (1/8" O.D. tube)
Cooling coil w or w/o valve (1/8" O.D. tube)
Process Thermocouple, Type J or K

External Accessories Available:

Vent Valve (1/8" O.D. tube)
2.5" (63.5mm) Dial Pressure Gauge - (Multiple ranges available)
Pressure transducers - range dependent on gauge
Inlet valves, (1/8" O.D. tube) either one or two on a shared connection
Catalyst charging valve 3/8" O.D. tube, 1/4" full bore opening
External thermocouple type J or K
Bottom port 3/16" (centered AE F437FB connection)

The following Engineering drawings are available upon request from Parker Autoclave Engineers for more detailed technical information.

Drawing Number 40A-8362 - Bench Top/Light Floor Motor Options (Air and DC motors)
Drawing Number 30B-0792 - Belt Drive Assembly (AC Motor)
Drawing Number 30A-9638 - Manual Screw Jack Assembly

Drawings				Drawing Title
316 Stainless Steel		Hastelloy® C-276		
100ml	300ml	100ml	300ml	
40C-0525	40A-8545	40C-0439	40A-8684	
	40A-8521		40A-8636	Bolted Closure Reactor
30A-9605	30A-9605	30B-0382	30A-0382	MAG 075 MagneDrive® Assembly
30A-9640	30A-9640	30B-0479	30A-0479	1/8" Valve Rack

Viton® and Kalrez® are registered trademarks of DuPont Dow Elastomers, Wilmington, DE
HASTELLOY® is a registered trademark of Haynes International Inc. Kokomo, IN

Note: Parker Autoclave Engineers, division of Snap-tite, Inc. reserves the right to substitute an equivalent material for trademarked material. Autoclave Engineers purchases substitute materials based on specification conformance, typically a widely accepted specification created by an industry standards organization.

Please refer to the following sections of the catalog for complimentary products and additional technical details. If your catalog is incomplete or out-of-date, feel free to register your name and download literature from Parker Autoclave Engineers web site. A good starting point is www.autoclaveengineers.com to reach the main page of Parker Autoclave Engineers reactor products.

Instrumentation - Details Parker Autoclave Engineers' full line of control options for temperature, pressure, and speed.

Agitation - Provides additional specifications on the MagneDrive® magnetic agitator and available impeller systems.

Pressure Vessels - Provides details on the Bolted Closure vessel assembly.

Stirred Reactor Selection Guide - Provides general information on all of Parker Autoclave Engineers' stirred reactors

Supporting Information

